

---

Subject: Omega Pro 18, Box Volume

Posted by [Kittikun](#) on Mon, 08 Jul 2002 23:41:34 GMT

[View Forum Message](#) <> [Reply to Message](#)

---

Hi Wayne, When I run Pialign, the recommend box volume is 3.3 cu.ft but from Pi web site, the box volume for stage 4 Pi-18 is the same as Thertter 4 Pi around 8 cu.ft. Therefore I am a little bit confuse. What is the suitable box volume for the Stage 4 Pi-18? Please. Best Regards, Kittikun