

---

Subject: 2 pi and PiAlign program  
Posted by [Paul](#) on Sun, 14 May 2006 13:09:17 GMT  
[View Forum Message](#) <> [Reply to Message](#)

---

the Qd value for the alpha 10 in the pialign data base that comes with the software is 2 but according tot he emminence website qts is .59 for the alpha 10 and the Qd should be  $1/.59 = 1.69$ . When I plug in 1.69 into pialign I get a slightly bigger enclosure that the one in the 2pi design. Why use a different QTS than what is provided in the emminence spec sheet? Thanks Paul

---